



Testing Equipment Manufacturer
For Paint, Coating & Ink

www.bevsinfo.com

BEVS Industrial Co.,Ltd. is the ISO approved laboratory testing instruments manufacturer in the field of coating, paint, resin&ink. Our quality control and testing instruments are made within the highest standards and combined with today's most advanced technology. BEVS will be the best partner with you and provide the whole solution of testing, inspection and laboratory instruments.



Viscosity

Cup Stand

Introduction:

BEVS Flow Cup Stand is used in concert with Ford Cup, DIN Cup, Afnor Cup and other similar cups, provides the level of platform for the cups. It is made up of three adjustable stainless steel poles and aluminum level bubble plate.



Technical Specification:

- ◆ Material: Poles—stainless steel or aluminum
Aluminum ring plate with level bubble

Order Information:

Order No.	Application Range
BEVS 1102/1	Ford Cup, Afnor Cup
BEVS 1102/2	DIN Cup