



Testing Equipment Manufacturer
For Paint, Coating & Ink

www.bevsinfo.com

BEVS Industrial Co.,Ltd. is the ISO approved laboratory testing instruments manufacturer in the field of coating, paint, resin&ink. Our quality control and testing instruments are made within the highest standards and combined with today's most advanced technology. BEVS will be the best partner with you and provide the whole solution of testing, inspection and laboratory instruments.



Drying

Drying Time Recorder


Introduction:

- ◆ BEVS Drying Time Recorder is used for the test of drying time or gel behavior of a variety of paints and coatings. It has the characteristics of intuition and good reproducibility, and can accurately evaluate drying time in all stages.
- ◆ Initiative design of touch screen control, customized user-defined start time and test time, make the operation more intuitive and simple.



Technical Specification:

- ◆ Dimension: 550×530×200 mm
- ◆ The number of tracks: 10
- ◆ Length of strips: 300 mm
- ◆ Driver: 5 sets of motors, 2 track / motor
- ◆ Drying Time:
 - Standard track: 6, 12, 24, 48 hours;
 - Customized track (track 1): Except the standard time, the user can customize the start and test time.
- ◆ Function reset during the testing: During the test, user can change the test time and re-start again.




Test Timing Table:

Testing Time(hour)						
Track No.	6	12	24		User Defined	Function Reset
1	√	√	√	√	√	√
2	√	√	√	√		√
3	√	√	√	√		√
4	√	√	√	√		√
5	√	√	√	√		√

Remark: "√": The test can be performed.

"User Defined": User can set the test time between 1 to 48 hours upon request, so that can conduct the personalized test.

"Function reset": During the test, user can change the test time and re-start again.



Standard:

ASTM D5895

Order Information:

BEVS1815 Drying Time Recorder