



Testing Equipment Manufacturer  
For Paint, Coating & Ink

www.bevsinfo.com

**BEVS Industrial Co., Ltd.** is the ISO approved laboratory testing instruments manufacturer in the field of coating, paint, resin & ink. Our quality control and testing instruments are made within the highest standards and combined with today's most advanced technology. BEVS will be the best partner with you and provide the whole solution of testing, inspection and laboratory instruments.



## Specific Gravity & Grinding

### Fineness of Grind Gauge

#### Introduction:

BEVS Fineness of Grind Gauge is manufactured according to the GB/T1724, ISO1524 standard.

#### Application:

Paint, plastic, pigment, printing ink, paper, ceramic, pharmaceutical, food and many other industries.



#### Technical Specification:

Material: Stainless steel

#### Order Information:

Item	Order No.	Effective Groove Size (L×W) (mm)	Measurement Range (μm)	Dimension (mm)	Division Value	Unit	Number of Grooves
Single-Channel Grind Gauge	1903/25	140×12.5	0-25	175×65×13	2.5	H & μm	1
	1903/50	140×12.5	0-50	175×65×13	5	H & μm	1
	1903/100	140×12.5	0-100	175×65×13	10	H & μm	1
Double-Channel Grind Gauge	1908/25	140×12.5	0-25	175×65×13	2.5	H & μm	2
	1908/50	140×12.5	0-50	175×65×13	5	H & μm	2
	1908/100	140×12.5	0-100	175×65×13	10	H & μm	2

#### Order Information of Accessories

BEVS1907	Scraper For Fineness of Grind Gauge
----------	-------------------------------------